

## SLIDE SCANNING INSTRUCTIONS

---

### INSTRUCTIONS

- ❖ Prepare samples
  - Must have a coverslip (#1.5 or 0.17mm)
  - Must be clean and fully dried
  - Circled sections to be imaged with a black marker: Not necessary if area to be scanned is very obvious (e.g. clean slide with single section and good H&E signal)
  - Clearly label slides
  - **Samples inappropriately prepared will NOT be scanned**
  
- ❖ Fill out “Slide Scanning Requisition” form (see attached)
  - One form per experimental set
  
- ❖ Drop off slides and form in room 778 before 10am
  - Slides should be placed in slide box
  
- ❖ Pick up
  - When samples are ready you will receive an email
  - Bring an external hard drive to download your data
  - **The Bioimaging facility is not responsible for backing up or storing data**
  - Don't forget to pick up your samples